

Assembling guide WM stretcher

1. Unpack the stretcher from the bag and unfold.



2. Assemble the aluminum pipes and put them into the longitudinal pockets.



3. Place the patient in the stretcher keeping the feet against the foot rest. Secure the feet with velcro straps.



4. Adjust the head rest to fit the patient.



5. Tie down the patient longitudinal with the belts between the patient's legs and secure above the shoulders.



6. First wrap the smaller flap then the bigger flap with airbag around the patient.

Secure the patient with the four cross belts.

Hold up the pipes close to the patient and tighten the cross belts.



7. Dress the patient with the neoprene hood.

8. Inflate the three airbags to stabilize the patient and for buoyancy in water.

9. Pull the spray hood over the face of the patient and upper airbag.

10. The stretcher is now ready for lifting either in vertical or horizontal position or in an angle of 45 degrees.

